

Veterinary ultrasound solution for Small and large animal pregnancy



WristScan V9

WristScan V9 is small in size, but big in functionality and features. With dust-proof and water-proof WristScan V9, you could make earlier, faster and more accurate diagnosis at the point of animal pregnancy at anytime, anywhere, no matter at animal hospital, even in or outdoor.



WristScan V9



Both specialty rod and probe are waterproof and dustproof, can be cleaned easily by flowing water.

- Designed for both small and large animal pregnancy — pigs, cows, sheep to dogs, cats and more
- High-quality images
- As easy to use, as it is easy to learn
- 340g (including the battery) light weight, small size and ergonomics designing
- More than 2.5 hours for continuous operation, ready anywhere, anytime with internal rechargeable battery and instant power-up
- Waterproof and dustproof probe, and can be easily cleaned by flowing water

Standard Configuration:

- One(1) main unit of V9 veterinary ultrasound system
- One(1) multi-frequency waterproof probe
- One(1) rechargeable Li Ion battery pack
- One(1) AC-adaptor (AC110V~240V, 50/60Hz)
- One(1) wrist bandage
- One(1) carrying case
- Two(2) years factory warranty (one year for probe)

Optional:

- Car Lighter
- Video Goggles
- Additional rechargeable Li Ion battery pack



Specifications:

- Monitor: 3.5 inch high resolution color TFT LCD, 640 x 480 pixels
- Probe: Multi-frequency mechanical sector probe
- Frequency: 2.5/3.5/5.0MHz
- Depth: Max 160mm
- Scan Angle: 78 degree
- Gray Shade: 256 gray levels
- Built-in Memory: up to 128 frames
- Common Measurement: Distance, Area, Volume, Heart Rate
- OB Measurement: Pig-HL Sheep/Goat-USD
Bovine-BL, SL, HL Cat-HD, BD

- AC-Adapter: Dog-BD, GSD, CRL, HD, BD
110V~240V, 50/60Hz
- Battery: Internal Li Ion battery pack
More than 2.5 hrs for continuous working
Fully charges batteries in 2hrs (75% charge in 1 hr)
- Weight: 340g including the battery
171mm×93mm×62mm
- Packing information: Dimension 0.32×0.38×0.12m
N.W.(2.4KG), G.W.(3KG)

Shenzhen Sunway Medical Device Co. Ltd.

4F General Bldg, 1st Industry Area, Guanlong village, Xili Town, Nanshan District, Shenzhen 518055, China

Tel: +86-755-86114956

Fax: +86-755-33203478

Web Site: www.chinasunway.com E-mail: export@chinasunway.com

Hotline: 400-666-2478

Specifications subject to change without notice

PN:2011012500EN01