



V8 HandScan



6.5 MHz Linear probe 5.0 MHz Micro-Convex 3.5 MHz Convex

Specifications :

- Screen Size: 5.5 inch, TFT color LCD
- Display Mode: B+, B+B, B+M, M, 4B
- Gray Shade: 256 levels
- Scan Depth: 70~240mm, adjustable
- Cine Loop: ≥ 400 frames
- Built-in Memory: 64 images
- Focus Position: 2 points adjustable
- Annotation: date, clock, name, Sex, ID, Age full screen editor.
- Battery Working Time: ≥ 3 hrs
- Battery Charging Time: ≥ 2hr
- Net Weight: 700g
- Package: metallic case in a carton box
- Image Reservation: up/down, left/right

- Image Processing: Coded color, GAMMA, Histogram And Image Smoother, N.G (near gain) adjustable, F.G (far gain) adjustable, O.G (overall gain) adjustable
- TCG:
- Language: English/Chinese switchable
- Measurement: Distance, Circumference, Area, Volume and Heart Rate
- OB Measurement: EQUINE---GSD, BOVINE---BL/SL/HL, SWINE---HL, SHEEP---USD, CAT---HD/BD, DOG---GSD/CRL/ HD /BD

Standard Configuration:

- One(1)Main unit
- One(1)6.5MHz Vet Rectal Linear Probe
- One(1)AC-Adaptor
- One(1)Rubber Protector
- One(1)Rechargeable Battery
- Two(2)shoulder and waist belts
- All above will be packed in a carrying case

Optional:

- Probes: R50/3.5MHz Convex Probe, L40/7.5MHz Linear Probe, R20/5.0MHz Micro-convex Probe
- Car Lighter
- Additional rechargeable battery
- Video Goggles
- Scanner Bracket



HandScan V8

Veterinary ultrasound scanner

HandScan V8 is a small, compact and battery powered UltraScanner for Veterinary applications anywhere in the field.

Shenzhen Sunway Medical Device Co. Ltd.

Hotline: 400-666-2478

4F General Bldg, 1st Industry Area, Guanlong village, Xili Town, Nanshan District, Shenzhen 518055, China

Tel: +86-755-86114956, Fax: +86-755-3203478

Web Site: www.chinasunway.com E-mail: export@chinasunway.com

specifications subject change without notice

PN:201008000EN