

CHAM CU-5000 (Black Type)

CE0123



Specifications

Specifications	
Power source	AC 220V, 50/60Hz
Power consumption	1200VA
Air compressor (Tank type)	250W, 2 way head max, air pressure 6Bar
Air compressor (No Tank type)	2.2Kgf/cm <sup>2</sup> ±0.7Kgf/cm <sup>2</sup>
Suction motor	Max. negative pressure 680mmHg Min. negative pressure 650mmHg
Ultraviolet sterilizer	253.7nm±5%
Dimension	1800(W) x 709(D) x 800(H)/ post pole : 1946(H)±10%
Weight	Approx 178kg±10%

Standard Set	Options
Spray	4 Microscope with light source & Arm
Mirror quick heating	1 Posr pole#2(Straight)
Ventilation	1 Monitor with arm & Base 1
Suction	2 Endoscopic system 1
100mm can	2 Doctor chair 1
83mm can	8 Patient chair 1
53mm can	3 Ear speculum stand 1
Waste bin	1 Irrigation 1
Instrument tray	1
Side can	1
Medicine bottle	11
Reflection light	1
Post pole#1(Curved)	1
UV Sterilization (Drawer type)	1

CHAM CU-5000 (White Type)

CE0123



Workstation

Patient chair

Endoscopic system

Microscope

Inhalation unit

Inspection equipment

Detail image



Two head suction motor



Suction bottle sliding holder



Ear&Nasal irrigation



Light on/off touch switch



Instrument tray&Heating tray



Key board holder



Adjustabel arm length



Centralizing main PCB



Main & Sub suction bottle



Side can



UV Sterilization (Drawer type)



Ear speculum stand