

XVS3 extreme

XVS5 extreme



세계일류상품



▲ XVS3 Extreme LCD cart type



세계일류상품



▲ XVS5 Extreme LCD cart type



▲ XVS3 Extreme



▲ CCD Gun type



▲ CCD Pen type



▲ XVS5 Extreme

Specification	Dimensions (WxDxH)(mm)	Weight	Light source	Camera	Resolution	Lamp life	Image transfer device
XVS3	314x303x110±10%	6.5kg±10%	LED 3.9V/13.5 one channel	CCD Gun type	768x494 pixels	60,000 hour	Optional
XVS5	314x303x160±10%	8.2kg±10%	LED 3.9V/13.5 two channel	CCD Pen type	768x494 pixels	60,000 hour	Optional

Standard set	CCD Camera	Camera adapter	Monitor	Moving cart	Main cable	Video cable	Light guide cable	Endoscope	Foot switch
XVS3	1ea	1ea	1ea	1ea	1ea	1ea	1ea	Otoscope(1ea)	1ea
XVS5	Optional	Optional	Optional	Optional	1ea	1ea	Optional	Optional	1ea

CONFIGURATIONS :

-Power source : AC 220V, 50/60Hz -Power consumption : 200VA