

GX-1

Electrical driven (Manual type)

Specifications

Power source	AC 220V, 50/60Hz
Power consumption	400VA
Rotation	340° Rotation by manual
Range of back-Rest movement	5°-95°
Up-Down stroke	200mm ±10%
Foot switch	Up / Down / Front / Rear / Reset (Automatic)
Dimension	600mm(ø)x1140mm(Min)x1350mm(Max)±10%
Weight	Approx 83kg
Armrest	Fixed

Standard set

Main body	1
Foot switch	1



Workstation

Patient chair

Endoscopic system

CH-200

Hydraulic driven

CE0123

Specifications

Power source	AC 220V, 50/60Hz
Power consumption	460VA
Rotation	340° Rotation by manual
Range of back-Rest movement	5°-95°
Up-Down stroke	230mm ±10%
Foot switch	Up / Down / Front / Rear / Reset (Automatic)
Dimension	600mm(ø)x1140mm(Min)x1350mm(Max)±10%
Weight	Approx 97kg
Armrest	360° Rotation by manual

Standard set

Main body	1
Foot switch	1



Microscope

Inhalation unit

Inspection equipment